



## Database Design and Description for CDR DICOM

Schick Technologies, Inc.  
30-00 47<sup>th</sup> Avenue  
Long Island City, NY 11101

(718) 937-5765  
(718) 937-5962 (fax)

Copyright © 2002 by Schick Technologies, Inc.  
All Rights Reserved

Many of the designations used by manufacturers and sellers to distinguish their products are claimed as trademarks. Where those designations appear in this document, and Schick Technologies, Inc. was aware of a trademark claim, the designations have been printed in caps or initial caps.

Part Number  
B1051050 Rev. –

August 15, 2002



Printed in the United States of America

This document was originally prepared in English

# Contents

<b>1. Introduction .....</b>	<b>1</b>
1.1. Intended Audience .....	1
1.2. CDR DICOM and SQL Databases .....	1
1.3. CDR DICOM Database Tables .....	1
1.3.1. Patient Table .....	2
1.3.2. Study Table .....	2
1.3.3. Series Table .....	2
1.3.4. Image Table .....	3
1.3.5. Volume Table .....	3
1.3.6. Workstations Table .....	4
1.3.7. Settings Table .....	4
<b>2. Step-by-Step Instructions.....</b>	<b>5</b>
2.1. Setting Up Sybase SQL Anywhere for CDR DICOM .....	5
2.2. Setting Up Microsoft SQL Server for CDR DICOM.....	6
<b>3. Sample SQL Scripts .....</b>	<b>7</b>
3.1. Sybase SQL Create Code .....	7
3.2. Microsoft SQL Server Create Code .....	10

# Notice

Schick Technologies shall not be liable for errors contained herein or consequential damages in connection with the furnishing, performance, or use of this document. We reserve the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes.

Copyright © 2002 Schick Technologies, Inc. All Rights Reserved

# 1. Introduction

## 1.1. Scope and Intended Audience

This document describes the content and format of database tables defined in CDR DICOM and the steps needed to incorporate these elements into Sybase and Microsoft SQL databases. The intended user of this document should be familiar with the design and function of SQL databases and their implementation in a client / server environment.

## 1.2. CDR DICOM and SQL Databases

The *CDRServer* application uses an SQL database to store information about patients and images and other CDR-related data objects. The database only contains information about images; the actual image file (in DICOM format) is stored separately on the server hard drive. This document describes the design of the database, including table and field descriptions as well as relationships between tables.

## 1.3. CDR DICOM Database Tables

The primary tables in the CDR Database are based on a simplified DICOM real world model as shown in **Figure 1**.




Figure 1. DICOM Real World Model

The corresponding tables in the CDR Database are **Patient**, **Study**, **Series** and **Image**. In addition, some supporting tables are **Volume**, **Workstations** and **Settings**. The **Volume** table stores information about specific storage volumes — hard drives, shared network drives, CD-ROMs, and other storage media. The **Workstations** table stores information about the workstations on the network that *CDRServer* recognizes. The **Settings** table is used to keep database version info and other useful information regarding the state of the database schema.

The fields for each of these tables are summarized on the following pages. In addition, **Figure 2** illustrates the one-to-many relationships among key fields in the tables.

### 1.3.1. Patient Table

Name	Type	Size	Description
PatientNumber	AutoNumber	4	Unique Identifier for Patient
NameDICOM	Text	64	DICOM formatted patient name LAST^FIRST
NameLast	Text	64	CDRPatient.LastName
NameFirst	Text	64	CDRPatient.FirstName
PatientID	Text	64	CDRPatient.IDNumber
PatientSex	Text	16	CDRPatient.Sex
BirthDate	Date/Time	8	CDRPatient.BirthDate
OtherID	Text	64	CDRPatient.OtherID
Comments	Memo	-	CDRPatient.Comments
DeleteFlag	Yes/No	1	Has patient been deleted?
EnteredDate	Date/Time	8	Date patient was entered
ModifiedDate	Date/Time	8	Date patient was last modified
PatientUID	Text	64	DICOM Unique ID for Patient

### 1.3.2. Study Table

Name	Type	Size	Description
StudyNumber	AutoNumber	4	Unique Identifier for Study
PatientNumber	Number (Long)	4	Link to Patient Table
StudyUID	Text	64	DICOM Unique ID of Study CDRStudy.UID
StudyDate	Date/Time	8	CDRStudy.Date
StudyID	Text	16	CDRStudy.StudyID
AccessionNumber	Text	16	CDRStudy.Accession
Description	Text	64	CDRStudy.Description
Comments	Memo	-	CDRStudy.Comments
ReferringPhysician	Text	64	CDRStudy.ReferringPhysician LAST^FIRST
ViewSetUID	Text	64	UID of Viewset Instance
ViewSetNumber	Number (Long)	4	Link to Image Table for Viewset entry
Exported	Yes/No	1	Has Study been exported?
Imported	Yes/No	1	Has Study been imported?
DeleteFlag	Yes/No	1	Has Study been deleted?
EnteredDate	Date/Time	8	Date Study was entered
ModifiedDate	Date/Time	8	Date study was last modified

### 1.3.3. Series Table

Name	Type	Size	Description
SeriesNumber	AutoNumber	4	Unique Identifier for Series
StudyNumber	Number (Long)	4	Link to Study Table
Modality	Text	16	CDRSeries.Modality
SeriesUID	Text	64	DICOM Unique ID of Series CDRSeries.UID
SeriesNumberDICOM	Text	12	Sequential number of Series in Study CDRSeries.Number
SeriesDate	Date/Time	8	CDRSeries.Date
PhysicianName	Text	64	CDRSeries.PhysicianName (formatted as LAST^FIRST)
OperatorName	Text	64	CDRSeries.OperatorName (formatted as LAST^FIRST)
SeriesDescription	Text	64	CDRSeries.Description
BodyPart	Text	16	CDRSeries.BodyPart
Exported	Yes/No	1	Has Series been exported?
Imported	Yes/No	1	Has Series been imported?
DeleteFlag	Yes/No	1	Has Series been deleted?
EnteredDate	Date/Time	8	Date Series was entered
ModifiedDate	Date/Time	8	Date Series was last modified

### 1.3.4. Image Table

Name	Type	Size	Description
ImageNumber	AutoNumber	4	Unique Identifier for Image
SeriesNumber	Number (Long)	4	Link to Series Table
ClassUID	Text	64	Class UID for DICOM Image
ImageUID	Text	64	DICOM Unique ID of Image CDRImage.UID
ImageNumberDICOM	Text	12	DICOM sequential number of image in series CDRImage.Number
ImageType	Text	255	CDRImage.ImageType
ImageFileName	Text	255	CDRImage.ImageFileName
ImageFileVolumeNumber	Number (Long)	4	Link to Volume Table for Image filename
ImageSource	Number (Long)	4	CDRImage.ImageSource
Width	Number (Long)	4	CDRImage.Width
Height	Number (Long)	4	CDRImage.Height
BitsPerPixel	Number (Long)	4	CDRImage.Depth
SensorType	Text	255	CDRImage.SensorType
SensorNumber	Text	16	CDRImage.SensorNumber
SensorMode	Text	255	CDRImage.SensorMode
MicronsPerPixel	Number (Long)	4	CDRImage.MicronsperPixel
MicronsPerPixelCalibrated	Number (Long)	4	CDRImage.MicronsPerPixelCal
ImageDate	Date/Time	8	CDRImage.ImageDate and ImageTime
AcquiredDate	Date/Time	8	CDRImage.AcqDate and AcqTime
CalibratedDate	Date/Time	8	CDRImage.CalibDate
Rotation	Number (Integer)	2	CDRImage.Rotation
HorizontalFlip	Yes/No	1	CDRImage.HorizontalFlip
Laterality	Text	16	CDRImage.Laterality
AnatomicRegion	Text	255	CDRImage.AnatomicRegion
AnatomicRegionModifier	Text	255	CDRImage.AnatomicRegionModifier
AnatomicStructure	Text	255	CDRImage.AnatomicStructure
AnatomicStructureModifier	Text	255	CDRImage.AnatomicStructureModifier
Comments	Memo	-	CDRImage.Comments
Exported	Yes/No	1	Has Image been exported?
Imported	Yes/No	1	Has Image been imported?
DeleteFlag	Yes/No	1	Has this Image been deleted?
EnteredDate	Date/Time	8	Date Image was entered
ModifiedDate	Date/Time	8	Date Image was last modified

### 1.3.5. Volume Table

Name	Type	Size	Description
VolumeNumber	Number (Long)	4	Unique ID of Volume
UNCPathName	Text	255	Drive Letter or UNC Share Name
CDRom	Yes/No	1	Is this volume a CD-ROM
Shared	Yes/No	1	Is this Volume shared
VolumeName	Text	255	Volume Label
VolumeSerial	Number (Long)	4	Volume serial number
FileSetID	Text	16	Name of File-Set Archive
FileSetUID	Text	64	Instance UID of File-Set Archive

### 1.3.6. Workstations Table

Name	Type	Size	Description
AETitle	Text	16	AE Title of workstation
IPAddress	Text	255	IP Address of workstation
Port	Text	255	Port used by workstation for listen
LocalAE	Text	16	AE Title used by server to initiate Association with workstation
MachineName	Text	255	Network machine name of workstation

### 1.3.7. Settings Table

Name	Type	Size	Description
Name	Text	255	Key Name of option
Type	Number (Long)	4	Data Type (0=String, 1=Number)
strValue	Text	255	String value of option
dwValue	Number (Long)	4	Numeric value of option




Figure 2. CDR Database Relationships



## 2. Step-by-Step Instructions

### 2.1. Setting Up Sybase SQL Anywhere for CDR DICOM

#### STEP 1

---

If you haven't done so already, please install the CDR DICOM Server application (CDRServer) before continuing with this procedure.

#### STEP 2

---

Using Sybase Central, create a new database.

#### STEP 3

---

Using Windows ODBC Administrator, create a new data source.

#### STEP 4

---

- A. Browse the CDR DICOM CD and copy the Sybase script found in the ODBC folder. (The script is also supplied in **Section 3.1.**)
- B. Start Interactive SQL and paste the script into the Command area. Click the **Execute** button.

#### STEP 5

---

- A. In CDR DICOM Server, click **Database > Configure**. At the Database tab, click on Data Sources to connect the server to the SQL database.
- B. Verify that your database connection is working correctly by importing patients, opening exams (after they've been added to the database), and making sure the database is able to connect to and function normally with registered workstations.

## 2.2. Setting Up Microsoft SQL Server for CDR DICOM

### STEP 1

---

If you haven't done so already, please install the CDR DICOM Server application (CDRServer) before continuing with this procedure.

### STEP 2

---

Using Enterprise Manager, create a new database.

### STEP 3

---

Using Windows ODBC Administrator, create a new data source.

### STEP 4

---

- A. Browse the CDR DICOM CD and copy the SQL Server script found in the ODBC folder. (The script is also supplied in **Section 3.2.**)
- B. Start SQL Query Analyzer and paste the script into the text area. Click the **Run** button (checkmark icon).

### STEP 5

---

- A. In CDR DICOM Server, click **Database > Configure**. At the Database tab, click on Data Sources to connect the server to the SQL database.
- B. Verify that your database connection is working correctly by importing patients, opening exams (after they've been added to the database), and making sure the database is able to connect to and function normally with registered workstations.

## 3. Sample SQL Scripts

### 3.1. Sybase SQL Create Code

Below is a sample SQL Script to create a CDR Database generated from Sybase SQL Modeler. It includes the creation of tables, default values, indexes, and foreign keys.

```
create table Settings
(
    Name                varchar(255)        not null,
    Type                smallint             ,
    strValue            varchar(255)        ,
    dwValue             integer              ,
    primary key (Name)
);

create table Workstations
(
    AETitle             varchar(16)          not null,
    IPAddress           varchar(255)        ,
    Port                varchar(255)        ,
    LocalAE             varchar(16)         ,
    MachineName         varchar(255)        ,
    primary key (AETitle)
);

create table Volume
(
    VolumeNumber        integer             not null
                        default AUTOINCREMENT,
    UNCPathName         varchar(255)        ,
    CDROM               BIT                not null,
                        default 0,
    Shared              BIT                not null,
                        default 0,
    VolumeName          varchar(255)        ,
    VolumeSerial        integer             ,
    FileSetID           varchar(16)         ,
    FileSetUID          varchar(64)         ,
    primary key (VolumeNumber)
);

create table Patient
(
    PatientNumber       integer             not null
                        default AUTOINCREMENT,
    NameDICOM           varchar(64)         not null,
    NameLast            varchar(64)         not null,
    NameFirst           varchar(64)         ,
    PatientID           varchar(64)         ,
    PatientSex          varchar(16)         ,
    BirthDate           DATETIME            ,
    OtherID             varchar(64)         ,
    Comments            long varchar        ,
    DeleteFlag          BIT                not null,
                        default 0,
    EnteredDate         DATETIME            not null,
                        default CURRENT_TIMESTAMP,
    ModifiedDate        DATETIME            not null,
                        default CURRENT_TIMESTAMP,
    PatientUID          varchar(64)         ,
    primary key (PatientNumber)
);

create index PatientNameID_IK on Patient (NameLast asc, NameFirst asc, PatientID asc);

create table Study
(
```

```

StudyNumber          integer          not null
    default AUTOINCREMENT,
PatientNumber        integer
    default AUTOINCREMENT,
StudyUID             varchar(64)      ,
StudyDate            DATETIME         ,
StudyID             varchar(16)       ,
AccessionNumber      varchar(16)      ,
Description          varchar(64)      ,
Comments            long varchar      ,
ReferringPhysician  varchar(64)      ,
ViewSetUID          varchar(64)      ,
ViewSetNumber       integer
    default AUTOINCREMENT,
Exported             BIT              not null
    default 0,
Imported            BIT              not null
    default 0,
DeleteFlag          BIT              not null
    default 0,
EnteredDate         DATETIME         not null
    default CURRENT_TIMESTAMP,
ModifiedDate        DATETIME         not null
    default CURRENT_TIMESTAMP,
primary key (StudyNumber)
);

```

```
create unique index StudyUID_IK on Study (StudyUID asc);
```

```
create table Series
```

```

(
SeriesNumber         integer          not null
    default AUTOINCREMENT,
StudyNumber         integer
    default AUTOINCREMENT,
Modality            varchar(16)       ,
SeriesUID           varchar(64)      ,
SeriesNumberDICOM   varchar(12)      ,
SeriesDate          DATETIME         ,
PhysicianName       varchar(64)      ,
OperatorName        varchar(64)      ,
SeriesDescription   varchar(64)      ,
BodyPart            varchar(16)      ,
Exported            BIT              not null
    default 0,
Imported            BIT              not null
    default 0,
DeleteFlag          BIT              not null
    default 0,
EnteredDate         DATETIME         not null
    default CURRENT_TIMESTAMP,
ModifiedDate        DATETIME         not null
    default CURRENT_TIMESTAMP,
primary key (SeriesNumber)
);

```

```
create unique index SeriesUID_IK on Series (SeriesUID asc);
```

```
create table Image
```

```

(
ImageNumber         integer          not null
    default AUTOINCREMENT,
SeriesNumber        integer
    default AUTOINCREMENT,
ClassUID            varchar(64)      ,
ImageUID            varchar(64)      ,
ImageNumberDICOM    varchar(12)      ,
ImageType           varchar(255)     ,
ImageFileName       varchar(255)     ,
ImageFileVolumeNumber
    integer
    default AUTOINCREMENT,
ImageSource         integer
    default 0,
Width               integer          ,
Height              integer          ,

```

```

BitsPerPixel          integer          ,
SensorType            varchar(255)      ,
SensorNumber          varchar(16)       ,
SensorMode            varchar(255)     ,
MicronsPerPixel       integer          ,
MicronsPerPixelCalibrated integer      ,
ImageDate             DATETIME         ,
AcquiredDate          DATETIME         ,
CalibratedDate        DATETIME         ,
Rotation              smallint         ,
HorizontalFlip        BIT              not null
    default 0,
Laterality            varchar(16)      ,
AnatomicRegion        varchar(255)     ,
AnatomicRegionModifier varchar(255)   ,
AnatomicStructure     varchar(255)     ,
AnatomicStructureModifier varchar(255) ,
Comments              long varchar     ,
Exported              BIT              not null
    default 0,
Imported              BIT              not null
    default 0,
DeleteFlag            BIT              not null
    default 0,
EnteredDate           DATETIME         not null
    default CURRENT_TIMESTAMP,
ModifiedDate          DATETIME         not null
    default CURRENT_TIMESTAMP,
primary key (ImageNumber)
);

create unique index ImageUID_IK on Image (ImageUID asc);

alter table Study
    add foreign key FK_STUDY_REF_102_PATIENT (PatientNumber)
        references Patient (PatientNumber);

alter table Series
    add foreign key FK_SERIES_REF_105_STUDY (StudyNumber)
        references Study (StudyNumber);

alter table Image
    add foreign key FK_IMAGE_REF_108_SERIES (SeriesNumber)
        references Series (SeriesNumber);

alter table Image
    add foreign key FK_IMAGE_REF_121_VOLUME (ImageFileVolumeNumber)
        references Volume (VolumeNumber);

```

## 3.2. Microsoft SQL Server Create Code

Below is a sample SQL Script to create a CDR Database compatible with Microsoft SQL Server. It includes the creation of tables, default values, indexes, and foreign keys.

```
create table Settings
(
    Name          varchar(255)      not null,
    Type          smallint          ,
    strValue      varchar(255)      ,
    dwValue       integer           ,
);

create table Workstations
(
    AETitle       varchar(16)       not null,
    IPAddress     varchar(255)      ,
    Port          varchar(255)      ,
    LocalAE       varchar(16)       ,
    MachineName   varchar(255)      ,
);

create table Volume
(
    VolumeNumber  integer           IDENTITY (1,1)    not null,
    UNCPathName   varchar(255)      ,
    CDRom        BIT               not null
    default 0,
    Shared        BIT               not null
    default 0,
    VolumeName    varchar(255)      ,
    VolumeSerial  integer           ,
    FileSetID     varchar(16)       ,
    FileSetUID    varchar(64)       ,
);

create table Patient
(
    PatientNumber integer           IDENTITY (1,1)    not null,
    NameDICOM     varchar(64)       not null,
    NameLast      varchar(64)       not null,
    NameFirst     varchar(64)       ,
    PatientID     varchar(64)       ,
    PatientSex    varchar(16)       ,
    BirthDate     DATETIME          ,
    OtherID       varchar(64)       ,
    Comments      text              ,
    DeleteFlag    BIT               not null
    default 0,
    EnteredDate   DATETIME          not null
    default CURRENT_TIMESTAMP,
    ModifiedDate  DATETIME          not null
    default CURRENT_TIMESTAMP,
    PatientUID    varchar(64)       ,
);

create index PatientNameID_IK on Patient (NameLast, NameFirst , PatientID );

create table Study
(
    StudyNumber   integer           IDENTITY (1,1)    not null,
    PatientNumber integer           ,
    StudyUID      varchar(64)       ,
    StudyDate     DATETIME          ,
    StudyID       varchar(16)       ,
    AccessionNumber varchar(16)     ,
    Description   varchar(64)       ,
    Comments      text              ,
    ReferringPhysician varchar(64)   ,
    ViewSetUID    varchar(64)       ,
    ViewSetNumber integer           ,
    Exported      BIT               not null
);
```

```

        default 0,
Imported          BIT          not null
        default 0,
DeleteFlag       BIT          not null
        default 0,
EnteredDate      DATETIME     not null
        default CURRENT_TIMESTAMP,
ModifiedDate     DATETIME     not null
        default CURRENT_TIMESTAMP,
);

create unique index StudyUID_IK on Study (StudyUID asc);

create table Series
(
    SeriesNumber      integer          IDENTITY (1,1)    not null,
    StudyNumber       integer          ,
    Modality          varchar(16)      ,
    SeriesUID         varchar(64)      ,
    SeriesNumberDICOM varchar(12)      ,
    SeriesDate        DATETIME         ,
    PhysicianName     varchar(64)      ,
    OperatorName      varchar(64)      ,
    SeriesDescription varchar(64)      ,
    BodyPart          varchar(16)      ,
    Exported          BIT              not null
        default 0,
    Imported          BIT              not null
        default 0,
    DeleteFlag       BIT              not null
        default 0,
    EnteredDate      DATETIME         not null
        default CURRENT_TIMESTAMP,
    ModifiedDate     DATETIME         not null
        default CURRENT_TIMESTAMP,
);

create unique index SeriesUID_IK on Series (SeriesUID asc);

create table Image
(
    ImageNumber       integer          IDENTITY (1,1)    not null,
    SeriesNumber      integer          ,
    ClassUID          varchar(64)      ,
    ImageUID          varchar(64)      ,
    ImageNumberDICOM varchar(12)      ,
    ImageType         varchar(255)     ,
    ImageFileName     varchar(255)     ,
    ImageFileVolumeNumber integer      ,
    ImageSource       integer          ,
    Width             integer          ,
    Height            integer          ,
    BitsPerPixel      integer          ,
    SensorType        varchar(255)     ,
    SensorNumber      varchar(16)      ,
    SensorMode        varchar(255)     ,
    MicronsPerPixel  integer          ,
    MicronsPerPixelCalibrated integer    ,
    ImageDate         DATETIME         ,
    AcquiredDate      DATETIME         ,
    CalibratedDate    DATETIME         ,
    Rotation          smallint        ,
    HorizontalFlip    BIT              not null
        default 0,
    Laterality        varchar(16)      ,
    AnatomicRegion    varchar(255)     ,
    AnatomicRegionModifier varchar(255) ,
    AnatomicStructure varchar(255)     ,
    AnatomicStructureModifier varchar(255) ,
    Comments          text             ,
    Exported          BIT              not null
        default 0,
    Imported          BIT              not null
        default 0,
    DeleteFlag       BIT              not null
);

```

```

        default 0,
EnteredDate          DATETIME          not null
        default CURRENT_TIMESTAMP,
ModifiedDate         DATETIME          not null
        default CURRENT_TIMESTAMP,
);

create unique index ImageUID_IK on Image (ImageUID asc);

ALTER TABLE [dbo].[Image] WITH NOCHECK ADD
    CONSTRAINT [PK_Image] PRIMARY KEY NONCLUSTERED
    (
        [ImageNumber]
    ) ON [PRIMARY] ,
    CONSTRAINT [IX_Image] UNIQUE NONCLUSTERED
    (
        [ImageUID]
    ) ON [PRIMARY]
GO

ALTER TABLE [dbo].[Patient] WITH NOCHECK ADD
    CONSTRAINT [PK_Patient] PRIMARY KEY NONCLUSTERED
    (
        [PatientNumber]
    ) ON [PRIMARY]
GO

ALTER TABLE [dbo].[Series] WITH NOCHECK ADD
    CONSTRAINT [PK_Series] PRIMARY KEY NONCLUSTERED
    (
        [SeriesNumber]
    ) ON [PRIMARY] ,
    CONSTRAINT [IX_Series] UNIQUE NONCLUSTERED
    (
        [SeriesUID]
    ) ON [PRIMARY]
GO

ALTER TABLE [dbo].[Study] WITH NOCHECK ADD
    CONSTRAINT [PK_Study] PRIMARY KEY NONCLUSTERED
    (
        [StudyNumber]
    ) ON [PRIMARY] ,
    CONSTRAINT [IX_Study] UNIQUE NONCLUSTERED
    (
        [StudyUID]
    ) ON [PRIMARY]
GO

ALTER TABLE [dbo].[Volume] WITH NOCHECK ADD
    CONSTRAINT [PK_Volume] PRIMARY KEY NONCLUSTERED
    (
        [VolumeNumber]
    ) ON [PRIMARY]
GO

ALTER TABLE [dbo].[Workstations] WITH NOCHECK ADD
    CONSTRAINT [PK_Workstations] PRIMARY KEY NONCLUSTERED
    (
        [AETitle]
    ) ON [PRIMARY]
GO

ALTER TABLE [dbo].[Settings] WITH NOCHECK ADD
    CONSTRAINT [PK_Settings] PRIMARY KEY NONCLUSTERED
    (
        [Name]
    ) ON [PRIMARY]
GO

GRANT SELECT , INSERT , UPDATE ON [dbo].[Image] TO [public]
GO

GRANT SELECT , INSERT , UPDATE ON [dbo].[Patient] TO [public]
GO

```



```

GRANT SELECT , INSERT , UPDATE ON [dbo].[Series] TO [public]
GO

GRANT SELECT , INSERT , UPDATE ON [dbo].[Study] TO [public]
GO

GRANT SELECT , INSERT , UPDATE ON [dbo].[Volume] TO [public]
GO

GRANT SELECT , INSERT , UPDATE, DELETE ON [dbo].[Workstations] TO [public]
GO

GRANT SELECT ON [dbo].[Settings] TO [public]
GO

ALTER TABLE [dbo].[Image] ADD
    CONSTRAINT [FK_Image_Series] FOREIGN KEY
    (
        [SeriesNumber]
    ) REFERENCES [dbo].[Series] (
        [SeriesNumber]
    )
GO

ALTER TABLE [dbo].[Image] ADD
    CONSTRAINT [FK_Image_Volume] FOREIGN KEY
    (
        [ImageFileVolumeNumber]
    ) REFERENCES [dbo].[Volume] (
        [VolumeNumber]
    )
GO

ALTER TABLE [dbo].[Series] ADD
    CONSTRAINT [FK_Series_Study] FOREIGN KEY
    (
        [StudyNumber]
    ) REFERENCES [dbo].[Study] (
        [StudyNumber]
    )
GO

ALTER TABLE [dbo].[Study] ADD
    CONSTRAINT [FK_Study_Patient] FOREIGN KEY
    (
        [PatientNumber]
    ) REFERENCES [dbo].[Patient] (
        [PatientNumber]
    )
GO

```